

Anti-RTL6 antibody

Catalog Number: 176055

Product name

Anti-RTL6 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to RTL6

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ldoc1l. aa range:141-190

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using LDOC1L Polyclonal Antibody