

Anti-SLC16A13 antibody



Catalog Number: 176068

Product name

Anti-SLC16A13 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to SLC16A13

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mot13. aa range:377-426

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 3T3 cells using MCT13 Polyclonal Antibody