

Anti-SLC16A14 antibody



Catalog Number: 176069

Product name

Anti-SLC16A14 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to SLC16A14

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mot14. aa range:131-180

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of HEPG2 cells using MCT14 Polyclonal Antibody diluted at 1: 1000