

Anti-C1orf35 antibody

Catalog Number: 176073

Product name

Anti-C1orf35 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to C1orf35

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mmtag2. aa range:10-59

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.;; Western blot analysis of 293T lysis using MMTAG2 antibody. Antibody was diluted at 1:1000