

Anti-MNT antibody

Catalog Number: 176075

Product name

Anti-MNT antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to MNT

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mnt. aa range:315-364

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using MNT Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).