

Anti-MRGPRX3 antibody



Catalog Number: 176087

Product name

Anti-MRGPRX3 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to MRGPRX3

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mrgx3. aa range:99-148

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of K562 cells using MRGX3 Polyclonal Antibody