

Anti-ABCC10 antibody

Catalog Number: 176089

Product name

Polyclonal

Anti-ABCC10 antibody

Ig Type**Specificity**

Rabbit IgG

Human, Monkey

Applications**Antibody description**

WB

Rabbit polyclonal antibody to ABCC10


Dilutions**Preparation**

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mrp7. aa range:835-884

Validations**Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Fig1:; Western Blot analysis of various cells using MRP7 Polyclonal Antibody diluted at 1: 1000

Clonality