

Anti-MRPL16 antibody



Catalog Number: 176091

Product name

Anti-MRPL16 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to MRPL16

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mrpl16. aa range:171-220

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using MRP-L16 Polyclonal Antibody diluted at 1: 2000