

Anti-MRPS34 antibody



Catalog Number: 176100

Product name

Anti-MRPS34 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to MRPS34

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mrps34. aa range:169-218

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using MRP-S34 Polyclonal Antibody diluted at 1: 1000