

# Anti-Atp1a1 antibody



Catalog Number: 176107

## Product name

Anti-Atp1a1 antibody

## Specificity

Rat

## Antibody description

Rabbit polyclonal antibody to Atp1a1

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from rat atp1 alpha1/na+k+ atpase1. aa range:15-64

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB, IF

## Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

## Validations




Fig1.; Western Blot analysis of HeLa cells using Na<sup>+</sup>/K<sup>+</sup>-ATPase α1 Polyclonal Antibody