

Anti-ANO7 antibody



Catalog Number: 176128

Product name

Anti-ANO7 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ANO7

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tm16g. aa range:881-930

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of 3T3 cells using NGEF Polyclonal Antibody diluted at 1: 1000