

Anti-UBAP2L antibody



Catalog Number: 176129

Product name

Anti-UBAP2L antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to UBAP2L

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ubap2l. aa range:561-610

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of NIH-3T3 K562 NIH-3T3 293T cells using NICE4 Polyclonal Antibody diluted at 1: 500