

Anti-NKX2-6 antibody

Catalog Number: 176130

Product name

Anti-NKX2-6 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to NKX2-6

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human nkx26. aa range:91-140

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using Nkx-2.6 Polyclonal Antibody diluted at 1: 500