

Anti-OR10D4P antibody

Catalog Number: 176142

Product name

Anti-OR10D4P antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR10D4P

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or10d4. aa range:161-210

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.: Western Blot analysis of various cells using Olfactory receptor 10D4 Polyclonal Antibody