

Anti-OR10Z1 antibody

Catalog Number: 176156

Product name

Anti-OR10Z1 antibody

Specificity

Human, Monkey

Antibody description

Rabbit polyclonal antibody to OR10Z1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or10z1. aa range:201-250

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.: Western Blot analysis of COS-7 cells using Olfactory receptor 10Z1 Polyclonal Antibody