

Anti-OR2A7 antibody



Catalog Number: 176166

Product name

Anti-OR2A7 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR2A7

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or2a4/2a7. aa range:231-280

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations



Fig1.: Western Blot analysis of various cells using Olfactory receptor 2A4/7 Polyclonal Antibody