

Anti-OR2AG2 antibody



Catalog Number: 176168

Product name

Anti-OR2AG2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR2AG2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or2ag1/2ag2. aa range:61-110

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of 293 cells using Olfactory receptor 2AG1/2 Polyclonal Antibody diluted at 1: 500