

Anti-OR2AJ1 antibody

Catalog Number: 176169

Product name

Anti-OR2AJ1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR2AJ1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or2aj1. aa range:201-250

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of Jurkat cells using Olfactory receptor 2AJ1 Polyclonal Antibody