

Anti-OR8K3 antibody

Catalog Number: 176235

Product name

Anti-OR8K3 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR8K3

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or8k3. aa range:236-285

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of LOVO cells using Olfactory receptor 8K3 Polyclonal Antibody