

Anti-OLFML2A antibody



Catalog Number: 176256

Product name

Anti-OLFML2A antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OLFML2A

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human olfml2a. aa range:311-360

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.: Western Blot analysis of RAW264.7 cells using Photomedin-1 Polyclonal Antibody