

Anti-PIP5K1C antibody

Catalog Number: 176258

Product name

Polyclonal

Anti-PIP5K1C antibody

Ig Type**Specificity**

Rabbit IgG

Human

Applications**Antibody description**

WB

Rabbit polyclonal antibody to PIP5K1C


Dilutions**Preparation**

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human pip5k1c. aa range:305-354

Validations**Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Fig1.; Western Blot analysis of various cells using PIPK I γ Polyclonal Antibody

Clonality