

Anti-PRDM12 antibody



Catalog Number: 176272

Product name

Anti-PRDM12 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to PRDM12

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human prdm12. aa range:191-240

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of various cells using PRDM12 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotec, MN, USA).