

Anti-RBSN antibody



Catalog Number: 176281

Product name

Anti-RBSN antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to RBSN

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human zfyve20. aa range:145-194

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using Rabenosyn-5 Polyclonal Antibody diluted at 1: 1000