

Anti-RPL10L antibody



Catalog Number: 176291

Product name

Anti-RPL10L antibody

Specificity

Human, Monkey

Antibody description

Rabbit polyclonal antibody to RPL10L

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rpl10l. aa range:151-200

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations



Fig1::; Western Blot analysis of various cells using Ribosomal Protein L10L Polyclonal Antibody diluted at 1: 1000