

# Anti-DENND2B antibody



Catalog Number: 176324

## Product name

Polyclonal

Anti-DENND2B antibody

## Ig Type

## Specificity

IgG

Human, Mouse

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to DENND2B

## Dilutions


## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human st5. aa range:411-460

## Validations

(kD)



## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Fig1.; Western Blot analysis of various cells using ST5 Polyclonal Antibody