

# Anti-RNF144A antibody



Catalog Number: 176366

## Product name

Anti-RNF144A antibody

## Ig Type

IgG

## Specificity

Human, Mouse

## Applications

WB

## Antibody description

Rabbit polyclonal antibody to RNF144A

## Dilutions

WB: 1:500-1:2000

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rnf144a. aa range:151-200

## Validations

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal



Fig1;; Western Blot analysis of various cells using UBCE7IP4 Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000