

Anti-USP50 antibody



Catalog Number: 176370

Product name

Anti-USP50 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to USP50

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human usp50. aa range:251-300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1:: Western Blot analysis of various cells using USP50 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000