

Anti-ATP6V1H antibody



Catalog Number: 176372

Product name

Anti-ATP6V1H antibody

Ig Type

IgG

Specificity

Human, Mouse

Applications

WB, IHC-P

Antibody description

Rabbit polyclonal antibody to ATP6V1H

Dilutions

WB: 1:500-1:2000

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human atp6v1h. aa range:341-390

IHC-P: 1:100-1:300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations



Fig1:: Western Blot analysis of various cells using V-ATPase H Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000