

Anti-ERCC5 antibody



Catalog Number: 176374

Product name

Anti-ERCC5 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ERCC5

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ercc5. aa range:131-180

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of various cells using XPG Polyclonal Antibody diluted at 1:2000. Secondary antibody (catalog# HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).