

Anti-ZADH2 antibody



Catalog Number: 176376

Product name

Anti-ZADH2 antibody

Ig Type

IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P

Antibody description

Rabbit polyclonal antibody to ZADH2

Dilutions

WB: 1:500-1:2000

Preparation

IHC-P: 1:100-1:300

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human zadh2. aa range:191-240

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations




Fig1:: Western Blot analysis of various cells using ZADH2 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000