

Anti-ZC3H13 antibody



Catalog Number: 176378

Product name

Anti-ZC3H13 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to ZC3H13

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human zc3h13. aa range:1571-1620

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using ZC3H13 Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000