

Anti-SLC39A9 antibody



Catalog Number: 176381

Product name

Anti-SLC39A9 antibody

Specificity

Human, Rat

Antibody description

Rabbit polyclonal antibody to SLC39A9

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human slc39a9. aa range:101-150

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using ZIP9 Polyclonal Antibody. Secondary antibody (catalog# HA1001) was diluted at 1:20000