

Anti-ZNF225 antibody



Catalog Number: 176383

Product name

Anti-ZNF225 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ZNF225

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human znf225. aa range:641-690

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations



Fig1.; Western Blot analysis of various cells using ZNF225 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).