

Anti-ZNF397 antibody



Catalog Number: 176387

Product name

Anti-ZNF397 antibody

Ig Type

IgG

Specificity

Human

Applications

WB, IHC-P, IF

Antibody description

Rabbit polyclonal antibody to ZNF397

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human znf397. aa range:10-59

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal




Fig1:: Western Blot analysis of various cells using ZNF397 Polyclonal Antibody. Secondary antibody □catalog#: HA1001) was diluted at 1:20000