

# Anti-ZNF682 antibody



Catalog Number: 176392

## Product name

Anti-ZNF682 antibody

## Ig Type

IgG

## Specificity

Human

## Applications

WB

## Antibody description

Rabbit polyclonal antibody to ZNF682

## Dilutions

WB: 1:500-1:2000

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human znf682. aa range:21-70

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Validations




Fig1.; Western Blot analysis of various cells using ZNF682 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).