

Anti-SLC30A9 antibody



Catalog Number: 176397

Product name

Polyclonal

Anti-SLC30A9 antibody

Ig Type

Specificity

IgG

Human, Mouse, Rat

Applications

Antibody description

WB

Rabbit polyclonal antibody to SLC30A9

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human slc30a9. aa range:181-230

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1.; Western Blot analysis of HepG2 cells using ZnT-9 Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000