

# Anti-APOL6 antibody



Catalog Number: 176400

## Product name

Anti-APOL6 antibody

## Ig Type

IgG

## Specificity

Human

## Applications

WB

## Antibody description

Rabbit polyclonal antibody to APOL6

## Dilutions

WB: 1:500-1:2000

## Preparation

This antigen of this antibody was synthesized peptide derived from apol6 . at aa range: 240-320

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Validations




Fig1:; Western Blot analysis of extracts from K562 cells, using ApoL6 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000