

Anti-APOL6 antibody



Catalog Number: 176400

Product name

Anti-APOL6 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to APOL6

Preparation

This antigen of this antibody was synthesized peptide derived from apol6 . at aa range: 240-320

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1:; Western Blot analysis of extracts from K562 cells, using ApoL6 Polyclonal Antibody.. Secondary antibody (catalog# HA1001) was diluted at 1:20000