

# Anti-SLC6A17 antibody



Catalog Number: 176406

## Product name

Anti-SLC6A17 antibody

## Ig Type

IgG

## Specificity

Human, Mouse, Rat

## Applications

WB

## Antibody description

Rabbit polyclonal antibody to SLC6A17


## Dilutions

WB: 1:500-1:2000

## Preparation

This antigen of this antibody was synthesized peptide derived from ntt4 . at aa range: 240-320

## Validations



## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

Fig1:; Western Blot analysis of extracts from rat stomach, using NTT4 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000