

# Anti-CCL4L1 antibody

Catalog Number: 176415

## Product name

Anti-CCL4L1 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to CCL4L1

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human ccl4. aa range:41-90

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1.; Western Blot analysis of HL-60, 22RV-1 cells using MIP-1 $\beta$  Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000