

Anti-NRG2 antibody



Catalog Number: 176416

Product name

WB, IHC-P

Anti-NRG2 antibody

Dilutions

Specificity

WB: 1:500-1:2000

Human, Mouse, Rat

IHC-P: 1:100-300

Antibody description

Rabbit polyclonal antibody to NRG2

Validations

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human nrg2. aa range:361-410

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications



Fig1:: Western Blot analysis of HT29 cells using Neuregulin-2 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at




Fig2:: Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:100