

Anti-RHAG antibody

Catalog Number: 176421

Product name

Anti-RHAG antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to RHAG

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the n-terminal region of human rhag. aa range:1-50

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-300

Validations




Fig1.; Western Blot analysis of K562 cells using CD241 Polyclonal Antibody.. Secondary antibody [catalog# HA1001) was diluted at 1:20000




Fig2.; Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100