

Anti-IGSF8 antibody

Catalog Number: 176437

Product name

Anti-IGSF8 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to IGSF8

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human igsf8. aa range:451-500

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of NIH-3T3 cells using CD316 Polyclonal Antibody.. Secondary antibody (catalog# HA1001) was diluted at 1:20000