

Anti-TNFRSF17 antibody



Catalog Number: 176465

General Information

Product name:

Anti-TNFRSF17 antibody

Immunogen:

Applications:

WB

Specificity:

Human

Antibody description:

Rabbit polyclonal antibody to TNFRSF17

Properties

Preparation:

This antigen of this antibody was synthesized peptide derived from the internal region of human cd269.

Formulation:

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality:

Polyclonal

Ig Type: IgG

Applications

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of L929 293T MCF7 HELA using CD269 antibody. Antibody was diluted at 1:2000. Secondary antibody (catalog# HA1001) was diluted at 1:20000