

# Anti-TDGF1P3 antibody

Catalog Number: 176467

## Product name

Anti-TDGF1P3 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to TDGF1P3

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human tdgf1p3. aa range:1-50

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1:; Western blot analysis of HepG2 293T AD293 HeLa lysis using TDGF1P3 antibody. Antibody was diluted at 1:1000. Secondary antibody (catalog# HA1001) was diluted at 1:20000