

# Anti-CMTM8 antibody



Catalog Number: 176468

## Product name

Polyclonal

Anti-CMTM8 antibody

## Ig Type

## Specificity

IgG

Human

## Applications

## Antibody description

WB, IHC-P

Rabbit polyclonal antibody to CMTM8

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human cmtm8. aa range:101-150

IHC-P: 1:100-1:300

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Validations

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality




Fig1:: Western blot analysis of hela lysis using CMTM8 antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000