

Catalog Number: 176483

---

## General Information

### Product name:

Anti-CCL24 antibody

### Immunogen:

### Applications:

IHC-P

### Specificity:

Human

### Antibody description:

Rabbit polyclonal antibody to CCL24

## Properties

### Preparation:

This antigen of this antibody was synthetic peptide from human protein at aa range: 61-110

### Formulation:

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

### Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### Clonality:

Polyclonal

### Ig Type: IgG

## Applications

IHC-P: 1:50-200

## Validations




Fig1.; Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200