

Anti-KCTD7 antibody



Catalog Number: 176501

Product name

Anti-KCTD7 antibody

Ig Type

IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

Rabbit polyclonal antibody to KCTD7

Dilutions

WB: 1:500-1:2000

Preparation

This antigen of this antibody was synthesized peptide derived from human kctd7. at aa range: 181-230

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations




Fig1:: Western Blot analysis of 1,mouse-brain 2,mouse-spleen cells using primary antibody diluted at 1:1000(4°C overnight). Secondary antibody: Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25°C, 1 hour)