

Anti-COL8A2 antibody



Catalog Number: 176503

Product name

Anti-COL8A2 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to COL8A2

Preparation

This antigen of this antibody was synthesized peptide derived from human col8a2. at aa range: 611-660

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1:; Western Blot analysis of 1,hela 2,3T3 cells using primary antibody diluted at 1:500(4°C overnight). Secondary antibody: Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25°C, 1 hour)