

# Anti-TMEM72 antibody



Catalog Number: 176507

## Product name

Anti-TMEM72 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to TMEM72

## Preparation

This antigen of this antibody was synthesized peptide derived from human tmem72 polyclonal

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1.; Western blot analysis of various lysate, antibody was diluted at 1000. Secondary antibody [catalog#: HA1001) was diluted at 1:20000